



HOVERMATT®

Mini (Half-Matt)


The HoverMatt® Mini (Half-Matt) delivers the same effortless lateral transfers and positioning as the standard HoverMatt® Air Transfer System, but in a shorter length. Specifically designed for use in the O.R. with specialty tables as well as in Labour and Delivery, the Mini is available in both reusable and single-patient use models. The reusable Mini with a double-coated finish is often preferred when excessive fluids make stain resistance a priority.



HoverMatt Mini (Half-Matt) Benefits:

- Air-assisted technology improves patient and caregiver safety and comfort during transfers
- Inflated mattress moves patient in stable position, reducing skin shear and bruising
- Half-length design is ideal for use with O.R. specialty tables
- Radiolucent and MRI compatible for artifact-free imaging
- Available in multiple styles and widths with a 544kg weight limit to accommodate most patients
- Variable speed air supply allows quieter and slower inflation options

HoverMatt Mini (Half-Matt)



For more information, contact your local HoverTech distributor, or Statina Healthcare International on 1300 791 404.

Product Specifications

Single Patient Use Mini (Half Matt)

Sewn Construction

Top: Non-woven polypropylene fibre

Bottom: Nylon fabric

Water soluble perimeter thread (120°)

Thread in sewn handles is NOT water soluble

Latex free

Available Sizes

86cm W x 114cm L - MHHT-HMSPU-34MI Box of 10

99cm W x 114cm L - MHHT-HMSPU-39MI Box of 10

Heat Sealed Mini (Half Matt)

Heat-Sealed Construction: RF welding in all seams

Material: Nylon twill

Latex free

Available Sizes

86cm W x 114cm L - MHHT-HMHS-34MI Each

Double Coated Mini (Half Matt)

Heat-Sealed Construction: RF welding in all seams

Material: Nylon twill

Latex free

Available Size:

86cm W x 114cm L - MHHT-HMDC-34MI Each

Compatible with HoverTech Air Supplies

MHHT-AIR2300 (Variable Speed - 230 V)

If a HoverMatt® Single Patient Use Mini (Half-Matt) is laundered, it will no longer inflate when connected to the air supply. The "DO NOT WASH" label at the foot-end of the product will dissolve and state "DISCARD IMMEDIATELY".

SHI221014-v2

Distributed By:

Imported into Australia by:
Statina Healthcare International
6/30 Leighton Place, Hornsby NSW
sales@statinainternational.com
Ph: 1300 791 404